

Growth hormone 1 Rat mAb, capture

Cat No.: SJ-AP2663

Applications

Reactive Species
Human

Positive Samples
N/A

Validated Applications
ELISA

Recommended Dilution
ELISA: 1/1000

Product Information

Antibody Type
Rat mAb

Immunogen
A synthesized peptide derived from GH1

Modification
Unmodified

Isotype
IgG1

Clone No.
N/A

Purification Method
N/A

Form
Liquid

Storage Buffer
Liquid in PBS containing 0.09% sodium azide as a preservative.

Storage Conditions
Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Host
Rat

Clonality
Monoclonal

Conjugation
Unconjugated

Concentration
N/A

Background Information

Function
N/A

UniProt ID
P01241

Organism
N/A

Recommended Name
N/A

Short Name
N/A

Alternative Names
GH1; GHN; Growth hormone 1; hGH-N; IGHD1B; Pituitary growth hormone; RNGHGP; Somatotropin

Gene ID
2688


Gene Name
GH1

Calculated MW
N/A

Observed MW
N/A

Subcellular Location
N/A

Validation Data



NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: support@sanjingbio.com

W: www.sanjingbio.com