

# Growth hormone 1 Rat mAb, unconjugated, detector



Cat No.: SJ-AP2664

## Applications

**Reactive Species**  
Human

**Positive Samples**  
N/A

**Validated Applications**  
ELISA

**Recommended Dilution**  
ELISA: 1/1000

## Product Information

**Antibody Type**  
Rat mAb

**Immunogen**  
A synthesized peptide derived from GH1

**Modification**  
Unmodified

**Isotype**  
IgG2a

**Clone No.**  
N/A

**Purification Method**  
N/A

**Form**  
Liquid

**Storage Buffer**  
Liquid in PBS containing 0.09% sodium azide as a preservative.

**Storage Conditions**  
Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

**Host**  
Rat

**Clonality**  
Monoclonal

**Conjugation**  
Unconjugated

**Concentration**  
N/A

## Background Information

**Function**  
N/A

**UniProt ID**  
P01241

**Organism**  
N/A

**Recommended Name**  
N/A

**Short Name**  
N/A

**Alternative Names**  
GH1; GHN; Growth hormone 1; hGH-N; IGHD1B; Pituitary growth hormone; RNGHGP; Somatotropin

**Gene ID**  
2688

**Gene Name**  
GH1


**Calculated MW**  
N/A

**Observed MW**  
N/A

**Subcellular Location**  
N/A

## Validation Data

---



**NOTICE:** For research use only. Not for therapeutic or diagnostic purposes.

## Contact

---

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: [support@sanjingbio.com](mailto:support@sanjingbio.com)

W: [www.sanjingbio.com](http://www.sanjingbio.com)