

# Syrian Hamster Anti-Mouse CD28 mAb (PE)



Cat No.: SJ-FC2825

## Applications

**Reactive Species**

Mouse

**Positive Samples**

N/A

**Validated Applications**

FC

**Recommended Dilution**

FC: 5  $\mu$ L per test for million cells or whole blood in 100  $\mu$ L volume.

## Product Information

**Antibody Type**

Syrian Hamster mAb

**Immunogen**

A synthesized peptide derived from Mouse CD28

**Modification**

Unmodified

**Isotype**

IgG

**Clone No.**

N/A

**Purification Method**

N/A

**Form**

Liquid

**Storage Buffer**

Liquid in PBS with 0.05% Proclin300, 1% BSA

**Storage Conditions**

Store at 2-8°C and protected from prolonged exposure to light. Avoid freeze/thaw cycles.

**Host**

Syrian Hamster

**Clonality**

Monoclonal

**Conjugation**

PE. Ex: 565nm. Em: 576nm.

**Concentration**

N/A

## Background Information

**Function**

N/A

**UniProt ID**

P31041

**Organism**

N/A

**Recommended Name**

N/A

**Short Name**

N/A

**Alternative Names**

T-cell-specific surface glycoprotein CD28; Cd28; CD28

**Gene ID**

12487

**Gene Name**

Mouse CD28

**Calculated MW**

N/A

**Observed MW**


N/A

**Subcellular Location**

N/A

## Validation Data

---



**NOTICE:** For research use only. Not for therapeutic or diagnostic purposes.

## Contact

---

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: [support@sanjingbio.com](mailto:support@sanjingbio.com)

W: [www.sanjingbio.com](http://www.sanjingbio.com)