

# HRP\* Goat Anti Mouse IgG(H+L)



Cat No.: SJ-SA4000

## Product Information

---

|              |                  |
|--------------|------------------|
| Host         | Reactive Species |
| Goat         | Mouse            |
| Applications | Conjugation      |
| WB,IHC,ELISA | HRP              |

## Detailed Information

---

### Recommended Dilution

WB 1:5000-20000;IHC 1:20-100;ELISA 1:5000-20000;

### Storage Buffer

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

### Storage Conditions

-15°C to -25°C/1 year(Do not lower than -25°C)

### Concentration

4 mg/ml

### Alternative Names

HRP\* Goat Anti-Mouse IgG;H+L

**NOTICE:** For research use only. Not for therapeutic or diagnostic purposes.

## Contact

---

For technical support and original validation data for this product please contact:  
T: 400-018-6916  
E: support@sanjingbio.com  
W: www.sanjingbio.com