

HRP* Goat Anti Mouse IgG(H+L)



Cat No.: SJ-SA4000

Product Information

Host Species

Goat

Reactive Species

Mouse

Validated Applications

WB,IHC,ELISA

Conjugation

HRP

Detailed Information

Recommended Dilution

WB 1:5000-20000;IHC 1:20-100;ELISA 1:5000-20000;

Storage Buffer

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage Conditions

-15°C to -25°C/1 year(Do not lower than -25°C)

Concentration

4 mg/ml

Alternative Names

HRP* Goat Anti-Mouse IgG;H+L

NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: support@sanjingbio.com

W: www.sanjingbio.com