

Rabbit Anti-Mouse IgG H&L (HRP)



Cat No.: SJ-SA6709

Product Information

Host	Reactive Species
Rabbit	Mouse
Applications	Conjugation
WB,ELISA	HRP

Detailed Information

Recommended Dilution

WB:1:5000-1:20000
ELISA:1:5000-1:20000

Storage Buffer

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Storage Conditions

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Concentration

N/A

Alternative Names

N/A

NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:
T: 400-018-6916
E: support@sanjingbio.com
W: www.sanjingbio.com