

Goat Anti-Rabbit IgG H&L (HRP) (Ready to Use)



Cat No.: SJ-SA6730

Product Information

Host	Reactive Species
Goat	Rabbit
Applications	Conjugation
IHC-P	HRP

Detailed Information

Recommended Dilution

Ready-to-use

Storage Buffer

liquid in PBS (without Mg and Ca) containing 150 mM □NaCl, 0.5% BSA, 0.02% sodium azide ,50% glycerol,and HRP-IgG polymer,PH7.4

Storage Conditions

Store at 4°C .

Concentration

N/A

Alternative Names

N/A

NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:

T: 400-018-6916

E: support@sanjingbio.com

W: www.sanjingbio.com