

Goat Anti-Chicken IgY H&L (Biotin)



Cat No.: SJ-SA6743

Product Information

Host	Reactive Species
Goat	Chicken
Applications	Conjugation
WB,IHC-F,IHC-P,ICC/IF,ELISA	Biotin

Detailed Information

Recommended Dilution

WB:1:5000-1:20000
IF:1:100-1:1000
ELISA:1:5000-1:20000

Storage Buffer

Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Storage Conditions

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Concentration

N/A

Alternative Names

N/A

NOTICE: For research use only. Not for therapeutic or diagnostic purposes.

Contact

For technical support and original validation data for this product please contact:
T: 400-018-6916
E: support@sanjingbio.com
W: www.sanjingbio.com